MUMMALIITICS

SPECIAL KIDS. SPECIAL NEEDS. SPECIAL SOLUTIONS.

PROFILE & PRODUCTS



- **9** Hope Island, QLD 4212
- 0476 806 917
- mummabilities@gmail.com
- www.mummabilities.com.au

PROUD DISTRIBUTOR OF



ABOUT MUMMABILITIES



I have been a special needs mum for more than 12 years now, so I would call myself a semi-pro with still so much to learn. Our beautiful son, Steele, was born with a very rare chromosome disorder, one that we only learnt about 4 days after his birth. I had a happy, healthy pregnancy, so this cyclone of lessons sure came as surprise. A surprise that I was somehow ready for. From the moment I kissed his warm lips, I knew I was his Mumma forever and that I, and his Dadda of course, were made to be his everything.

And so, our journey began. From doctors to specialists to surgeries, Steele has been (and still is) our greatest teacher.

At some stage, we all ponder our purpose in life. I believe a big part of mine is to reach out to those travelling a similar journey to ours and share what I have learnt, starting with resources and tools that have not only added a tremendous amount of value to Steele's progression, but made our family's life one little step easier each time.

Being a full-time working mum too, time was never on my side, but 2020 brought with it a new slate. I held my breath, jumped....and Mummabilities was born.

I am Steele's "Mumma" and I want to share with you my "abilities" – the lessons I have learnt with the means and skills I have adopted to help us along our journey.

Tarryn XXX

EXCLUSIVE DISTRIBUTOR

EASY UP



This is the product that made me realise I needed to put my ideas into action. My son has been wearing AFOs since the age of 2. We were introduced to these AFO friendly shoes over 6 years ago and haven't found any other product remotely similar in their functionality, comfort and style. Mummabilities is the exclusive distributor of EASY UP® shoes in Australia. I wanted to introduce these incredible hand-made Italian shoes to the Australian market and help make managing AFOs a little easier.

THE DESIGN



Handmade in Italy, the EASY UP® collection is designed for independent use by AFO, DAFO, splint users, or those with reduced or lost mobility.

Thanks to its "trap door like" rear opening, the shoe can completely open, allowing the foot (with AFO) to easily slide into place.

The innovative EASY UP BOA® lacing system closes the "trap door" and adjusts the shoe to fit the foot perfectly. This means that there are no pressure points on the instep that tend to be common with laces and velcro. The simplicity of the EASY UP BOA® lacing system means that the shoe can be placed on and off with one hand, or easily by the carer.

CUSTOMISE YOUR WAY



EASY UP® shoes are custom made to your selection. All you have to do is:

- 1. Decide on a suitable dial position (side or back)
- 2. Select your preference in model/style (Children's, Men's or Ladies)
- 3. Select your preference in lining, sole and heel
- 4. Verify your size
- 5. Place your order online via Paypal, request an invoice or email

CUSTOM MADE



PRODUCT RANGE

The EASY UP® product range has a style to suit everyone. There are 20 styles in the Children's collection and 5 styles in both the Women's and Men's collections.

There are 2 position options for the dial of the EASY UP BOA® lacing system. Some EASY UP® models have a rear dial and some have a side dial, which would better suit wheelchair users or those with restricted movement.



SIDE & BACK DIAL MODELS

CHILDREN'S SIDE DIAL MODELS - Air, Buddy, Grow, Kid, Neo, Olly (sandal), Quattro, Quattro (sandal), Quilt, Rain, Silver, Sock, Snow and Wing.

CHILDREN'S BACK DIAL MODELS - Explorer, Ice, Laser, Rover, Star, Trekker (winter).



ADULT SIDE DIAL MODELS - Mens: WMS and WMB. Ladies: WWS and WWB.

ADULT BACK DIAL MODELS - Mens: SM, SMB and FM.

Ladies: SW, SWB and FW.



LINING, SOLES & HEELS

When selecting the options best suited for you or the person you care for, the following can be considered:

INNER LINING - Plain (full) or perforated lining for better breathability.

SOLES - Standard (firm), Flexible or Special walker sole.

HEELS - Flat heel (greater stability) or rolling heel (normal gait/DAFOs).

CHILDREN'S RANGE

AIR COLLECTION







BUDDY COLLECTION



















CHILDREN'S RANGE

BUDDY COLLECTION

























CHILDREN'S RANGE

BUDDY COLLECTION

























CHILDREN'S RANGE

BUDDY COLLECTION



STAR COLLECTION

















CHILDREN'S RANGE

KID COLLECTION





















CHILDREN'S RANGE

LASER COLLECTION









ICE COLLECTION









CHILDREN'S RANGE

SILVER COLLECTION













NEO COLLECTION





CHILDREN'S RANGE

QUATTRO COLLECTION





QUILT COLLECTION







RAIN COLLECTION





CHILDREN'S RANGE

SNOW COLLECTION









SOCK COLLECTION







CHILDREN'S RANGE

WING COLLECTION









GROW COLLECTION





QUATTRO SANDAL COLLECTION





CHILDREN'S RANGE

OLLY SANDAL COLLECTION







TREKKER COLLECTION





ROVER COLLECTION





EXPLORER COLLECTION





MEN'S RANGE

WALK COLLECTION





WALK BOOT COLLECTION





SMART COLLECTION





SMART BOOT COLLECTION





MEN'S RANGE

FREE COLLECTION

FM 01



FM 02



FM 03



FM 04



FM 05



LADIES RANGE WALK COLLECTION





WALK BOOT COLLECTION





SMART COLLECTION





SMART BOOT COLLECTION





LADIES RANGE

FREE COLLECTION











SIZE CHART



VERIFY YOUR SIZE

Sizing

Please use the chart below as a reference when verifying the correct size:

A: OVERALL LENGTH

B:TWIST TOES PERIMETER (circumference)


B1: TWIST TOES WIDTH


C: METATARSAL PERIMETER / ANATOMIC BALL (circumference)

D: HIGH INSTEP / HEEL PERIMETER (circumference)

REFERENCE GUIDE

The following measurements are intended as a guide to assist you in determining the best fit in the Easy Up Kids range. They are based on production measurements. Tolerance \pm is allowed. Should your measurements largely deviate from one or more of the below references, please request assistance from your professional orthotist. A basic "rule" to follow (for Kids and Adults range) is to measure the length of your foot/feet with AFO on and reference this to A (Overall length).





SIZE CHART



All measurements must be taken while wearing AFO or splints. Please email mummabilities@gmail.com with any queries.

6177					
SIZE	A	В	С	D	B1
20	13,30 cm	16,70 cm	18,00 cm	22,20 cm	5.92 cm
21	13,95 cm	17,15 cm	18,50 cm	22,85 cm	6,10 cm
22	14,60 cm	17,60 cm	19,00 cm	23,50 cm	6,28 cm
23	15,30 cm	18,05 cm	19,50 cm	24,15 cm	6,46 cm
24	15,95 cm	18,50 cm	20,00 cm	24,75 cm	6,64cm
25	16,60 cm	18,95 cm	20,50 cm	25,40 cm	6,82 cm
26	17,30 cm	19,40 cm	21,00 cm	26,10 cm	7,00 cm
27	17,95 cm	19,85 cm	21,50 cm	26,65 cm	7,18cm
28	18,60 cm	20,30 cm	22,00 cm	27,30 cm	7,36cm
29	19,60 cm	20,80 cm	22,70 cm	28,05 cm	7,54cm
30	19,95 cm	21,25 cm	23,30 cm	28,70 cm	7,72 cm
31	20,60 cm	21,70 cm	23,70 cm	29,35 cm	7,90 cm
32	21,30 cm	22,15 cm	24,20 cm	30,00 cm	8,08 cm
33	21,95 cm	22,60 cm	24,70 cm	30,65 cm	8,26cm
34	22,55 cm	23,10 cm	25,50 cm	31,50cm	8,44 cm
35	23,10 cm	23,55 cm	26,00 cm	32,15 cm	8,62 cm
36	23,90 cm	24,00 cm	26,50 cm	32,85 cm	8,80 cm
37	24,55 cm	24,45 cm	27,00 cm	33,50cm	8,98 cm
38	25,25 cm	24,90 cm	27,60 cm	34,20cm	9,16cm
39	25,90 cm	25,35 cm	28,10 cm	34,90 cm	9,34cm
40	26,55 cm	25,80 cm	28,60 cm	35,55 cm	9,52 cm
41	27,25 cm	26,25 cm	29,20 cm	36,20cm	9,70 cm
42	27,90 cm	26,70 cm	29,70 cm	36,90 cm	9,88 cm
43	28,55 cm	27,15 cm	30,20 cm	37,55 cm	10,06 cm
44	29,25 cm	27,60 cm	30,70 cm	38,25 cm	10,24cm
45	29,90 cm	28,06 cm	31,30 cm	38,90 cm	10,42 cm
46	30,55 cm	28,50 cm	31,80 cm	39,55 cm	10,68 cm
47	31,25 cm	29,00 cm	32,30 cm	40,25 cm	10,86 cm

LINING, HEEL & SOLE

LINING OPTIONS



PLAIN



PERFORATED

HEEL OPTIONS



FLAT

The flat heel increases stability for those who do not walk independently and / or for those who wear non-dynamic splints.



ROLLING

The rolling heel helps those who walk independently with a normal gait (stance and swing phase) and / or those who wear dynamic splints / DAFO.

SOLE OPTIONS



STANDARD (FIRM)



FLEXI

LINING, HEEL & SOLE

SOLE OPTIONS continued



SPECIAL WALKER

Designed to accommodate walking aids/equipment. Please contact me to enquire about this option.

HOW TO ORDER

3 OPTIONS

- 1. Order online via Paypal at <u>www.mummabilities.com.au</u>
- 2. Request an invoice via https://www.mummabilities.com.au/ndis-form
- 3. Email your order to mummabilities@gmail.com

If you have any further questions, please contact me via email at mummabilities@gmail.com or call 0476 806 917.

MUMMALIITICS

SPECIAL KIDS. SPECIAL NEEDS. SPECIAL SOLUTIONS.



www.mummabilities.com.au



mummabilities@gmail.com



0476 806 917

PROUD DISTRIBUTOR OF

